

# FIND YOUR WAVELENGTH

UV

VISIBLE

SWIR

MWIR / LWIR

C-BLUE One  
CMOS  
400 – 1000 nm




OCAM<sup>2</sup>  
EMCCD  
400 – 1000 nm



C-BLUE 2 (Q4 2022)  
High Resolution CMOS  
400 – 1000 nm



C-RED One  
MCT e-APD  
800 – 2500 nm



C-RED 2  
InGaAs  
900 – 1700 nm



C-RED 2 Lite  
InGaAs  
900 – 1700 nm



C-RED 2 ER 2.2µm  
Extended InGaAs  
1300 – 2200 nm



C-RED 2 ER 1.9µm  
Extended InGaAs  
1100 – 1900 nm



C-RED 3  
InGaAs  
900 – 1700 nm



UV = Ultraviolet  
VIS = Visible  
SWIR = Short Wave InfraRed  
MWIR = Middle Wave InfraRed  
LWIR = Long Wave InfraRed